

The advanced solution for faster, effective wound care
Portable | Efficient | Easy-to-Us



VIVEREX

A negative pressure
wound therapy system

About Our NPWT Devices

Our **Single-Use Negative Pressure Wound Therapy** (sNPWT) devices are meticulously designed to accelerate the healing process for a broad spectrum of wounds. We offer two specialized variants: **VivereX™** and **VivereX+™**, both engineered to provide effective treatment across diverse wound types.

Two Variants

VivereX™ - Offers treatment durations of **10 days**



VivereX is a canister-free sNPWT system that combines quality and ease of use with affordability, providing efficient wound care solutions.

VivereX+™ - Offers treatment duration of **21 days**



VivereX+™ is the practical alternative for sNPWT providing more value through extended care. The treatment duration of 21 days is consistent with healing time for most lower extremity wounds and back.

Ideal for

The **VivereX™** & **VivereX+™** Single-use canister-less NPWT Systems are suitable for use in hospitals and clinics as well as a home care setting. They are designed to remove small to moderate volumes of exudate from below listed indications:

- **Acute wound**
- **Chronic wound with small or medium secretion**
- **Skin graft or flaps and donor site**
- **Pressure ulcers**
- **Diabetics ulcers**
- **Subacute dehisced wounds**
- **Surgical clean closed incision**
- **Trauma wounds**
- **Burn injury grade I and II**



Benefits of Single-Use NPWT

Single-use NPWT offers several benefits compared to traditional NPWT systems, making it a convenient and effective option for wound care.




Fewer Hospital Visits

Can be managed in outpatient settings, reducing hospital stays and healthcare costs.



Lower Upfront Cost

Less expensive than traditional NPWT units, making it a cost-effective option.



Lightweight & Compact

Easier for patients to carry and use in outpatient or home settings.




Better Mobility & Comfort

Smaller, quieter devices improve patient comfort and mobility.

How VivereX™ Works

Our sNPWT device uses **controlled negative pressure** to draw fluid away from the wound, creating an optimal healing environment. By applying consistent suction, the device promotes blood flow, reduces edema, and helps remove exudate from the wound area.



The VivereX(+) features a battery compartment, belt clip for easy attachment, and a tube connection for seamless therapy. A centrally placed power button ensures simple operation.

For optimal performance, the VivereX(+) features intelligent indicators in a compact design, ensuring reliable and effective wound care while remaining user-friendly for both patients and healthcare professionals.

Why VivereX™ & VivereX+™

1.

Canister-free system

Lighter, more compact & easier to hide.
Reduced risk of contamination.

2.

Silent

With a sound level of <50db
and no vibrations.

3.

Possibility of 21 days treatment

VivereX™ will cover the need of treatment for many wounds but with VivereX+,™ we introduce the possibility of 21 days of treatment.

4.

Low battery consumption

There is no need to replace batteries during treatment.

5.

Non-magnetic pump

The non-magnetic pump operates safely without affecting nearby pacemakers.

6.

Competitive Pricing Strategy

We are continuously aiming to provide the best possible solution for a competitive price.



Discover the Future of Wound Care.

Request a Demo Today!

E-mail: info@oscartech-uk.com


Phone: +44 2036338413

Website: www.oscartech-uk.com

Accelerate Wound Healing with Our Advanced NPWT Dressings

Our NPWT Dressings are designed to optimize wound care by enhancing tissue re generation, reducing infection risk, and minimizing healing time.

W e have 10 sizes available where 2 are flex dressings for difficult anatomic positions like elbow, knee and heel.



Dressing Sizes	
1	10 x 20 CM
2	10 x 30 CM
3	10 x 40 CM
4	16 x 16 CM
5	16 x 21 CM
6	16 x 31 CM
7	21 x 21 CM
8	26 x 26 CM

1. Controlled Negative Pressure

2. Easy to Apply, Secure Seal

3. High Absorbency

4. Even Negative Pressure

VivaFoam™

VivaFoam™ is a premium filler designed specifically to enhance the wound healing process when used in conjunction with the VivereX(+)[™] system. This innovative foam effectively removes excess exudates, infectious materials, and tissue debris from the wound site, creating an optimal healing environment.

It is recommended to use VivaFoam together with VivereX(+) when there is a wound with a depth of more than 5 mm.



VivaFoam Sizes

1	10 x 7.5 x 3.2 CM
2	15 x 10 x 1.5 CM
3	20 x 10 x 3.2 CM
4	25 x 15 x 1.5 CM
5	25 x 15 x 3.2 CM

Key Features of VivaFoam™:

- **Hydrophobic Design:** Repels moisture and fluids to maintain an ideal wound environment.
- **Open Cell Structure:** Promotes effective fluid management and maintains consistent negative pressure for enhanced therapy outcomes.
- **Soft and Conformable:** Conforms to the wound bed for maximum contact and coverage.
- **High transfer Capacity:** Transfer exudate effectively to minimize the risk of maceration.
- **Gentle on wounds:** Designed for patient comfort and minimal trauma during dressing changes.

Order Information

VivereX™				VivereX™ Dressings			
Part-nr	Item	Unit	Qty/unit	Part-nr	Item	Unit	Qty/Unit
PDI600121	VivereX with 2pcs 10x20 dressing	Single	1	PDI605221	Dressing 10x20	Case	5
PDI600122	VivereX with 2pcs 10x30 dressing	Single	1	PDI605222	Dressing 10x30	Case	5
PDI600123	VivereX with 2pcs 10x40 dressing	Single	1	PDI605223	Dressing 10x40	Case	5
PDI600124	VivereX with 2pcs 16x16 dressing	Single	1	PDI605224	Dressing 16x16	Case	5
PDI600125	VivereX with 2pcs 16x21 dressing	Single	1	PDI605225	Dressing 16x21	Case	5
PDI600126	VivereX with 2pcs 16x31 dressing	Single	1	PDI605226	Dressing 16x31	Case	5
PDI600127	VivereX with 2pcs 21x21 dressing	Single	1	PDI605227	Dressing 21x21	Case	5
PDI600128	VivereX with 2pcs 26x26 dressing	Single	1	PDI605228	Dressing 26x26	Case	5
VivereX+™				VivaFoam™			
PDI600421	VivereX+ with 2pcs 10x20 dressing	Single	1	PDI636301	VivaFoam 10x7.5x3.2	Case	36
PDI600422	VivereX+ with 2pcs 10x30 dressing	Single	1	PDI672305	VivaFoam 15x10x1.5	Case	72
PDI600423	VivereX+ with 2pcs 10x40 dressing	Single	1	PDI624302	VivaFoam 20x10x3.2	Case	24
PDI600424	VivereX+ with 2pcs 16x16 dressing	Single	1	PDI648303	VivaFoam 25x15x1.5	Case	48
PDI600425	VivereX+ with 2pcs 16x21 dressing	Single	1	PDI624304	VivaFoam25x15x3.2	Case	24
PDI600426	VivereX+ with 2pcs 16x31 dressing	Single	1	PDI600301	Vivafoam Filler 10 x 7.5 x 3.2	Single	1
PDI600427	VivereX+ with 2pcs 21x21 dressing	Single	1	PDI600302	Vivafoam Filler 20 x 10 x 3.2	Single	1
PDI600428	VivereX+ with 2pcs 26x26 dressing	Single	1	PDI600303	Vivafoam Filler 25 x 15 x 1.5	Single	1
				PDI600304	Vivafoam Filler 25 x 15 x 3.2	Single	1
				PDI600305	Vivafoam Filler 15 x 10 x 1.5	Single	1